

WEIGH.
WEIGH.
WEIGH.
WEIGH.



MEASURE.
MEAS **DETECTO** RE.
MEASURE.
MEASURE.



STORE.
MEAS **DETECTO** RE.
STORE.
STORE.



DISPOSE.
DISPOSE.
DISPOSE.
DISPOSE.






Weigh. Measure. Store. Dispose.

DETECTO is a leading manufacturer of precision medical and industrial weighing solutions, trusted by healthcare providers, laboratories, and professionals worldwide.

With more than a century of manufacturing expertise, DETECTO designs and builds scales, balances, indicators, medical carts, and more that deliver accuracy, durability, and ease of use.

From NTEP-certified medical scales to robust industrial products, DETECTO combines in-house engineering, rigorous quality standards, and thoughtful design to support reliable measurement in critical environments every day.



01

Physician Scales
solo®, apex®, icon®,
Weigh Beam

05

Medical Carts
Rescue, MobileCare,
Pediatric Scale

08

Stretcher Scales
8500

09

Wheelchair Scales
7150-AC, 6550, 6560,
BRW1000, 6400, 7550

15

Bariatric Scales
6855, 6856, 6868, 6800

19

Chair Scales
6475, 6880

21

Pediatric Scales
6745, 8432-CH, 8440, MB130

25

Stadiometers
DLM, PHR, SONARIS

28

In-Floor Dialysis Scales
Solace

29

Home Healthcare Scales
APEX-RI, DR660, GP-600-MV1,
DR400C/DR550C, SlimTALK, SlimPRO

35

Wet Diaper/Chux Scales
WPS12UT

36

Precision Digital
Balances
420

37

In-Bed Scales
IB

38

Patient Lift Scales
PL

39

Waste Receptacles
Step-On Cans, Waste Mate


41

Hanging Scales
T3530

42

Glove Box Holders
GH

solo® Eye-Level Digital Physician Scales with Inline Mechanical Height Rods



Inline mechanical height rod is situated to the side of the LCD, so it never interferes with viewing measurements.




Capacity	550 lb x 0.2 lb / 250 kg x 0.1 kg
Platform Size	14 in x 15 in x 2.4 in / 36 cm x 38 cm x 6 cm
Power	6 AA batteries (included) or AC adapter (included)
Height Rod Range	30–84.5 in / 76–214.5 cm
Other Features	Calculates weight, height, and BMI Integral wheels ease transport

- Most economical, high-quality digital physician scale with height rod
- Low-profile platform and sturdy column for stability
- Eye-level LCD display with 1-in-high (2.5-cm) readouts
- RS232 serial port for EMR/EHR data output
- Indicator features 8 simple buttons that are easily identifiable in any language

DETECTO®



apex® Series Digital Physician Scales with Mechanical or Sonar Height Rods



The apex® features a high-tech, smart-phone-style indicator for easy utilization with 6 simple buttons and a bold LCD display for large readouts.



- Extra wide and flat platform
- Choice of sonar touchless or inline mechanical height rod
- Touchless sonar height rod automatically locks on patient's height and measures nearly any patient size
- RS232 serial port and micro USB-B port for data transfer
- Available models with Wi-Fi/Bluetooth for wireless EMR/EHR and Welch Allyn® device connectivity (CVSM/CSM, LXI)



Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



600 lb x 0.1 lb /
270 kg x 0.05 kg

Capacity

17 in x 17 in x 2.75 in /
43.2 cm x 43.2 cm x 7 cm

Platform Size

12 AA batteries (not included)
or AC adapter (included)

Power

24–86 in x 0.125 in /
60–218 cm x 0.1 cm

**Height Rod
Range**

Mother/Baby function factors
the weight of infants and
toddlers held by an adult

Other Features

Integral wheels ease transport

Optional handrail

OneWeigh units locking
into pounds or kilograms

Calculates weight, height,
and BMI



icon® Eye-Level Digital Physician Scales with Touchless Sonar Height Rods



Capacity	800 lb x 0.1 lb / 400 kg x 0.05 kg
Platform Size	17 in x 17 in x 1.5 in / 43.2 cm x 43.2 cm x 3.8 cm
Power	12 AA batteries (not included) or AC adapter (included)
Height Rod Range	24–86 in x 0.2 in / 60–218 cm x 0.5 cm
Other Features	Displays weight, height, and BMI onscreen simultaneously Multi-color touchscreen LCD OneWeigh units locking into pounds or kilograms

DETECTO



Multi-color touchscreen enables intuitive entry of 8-digit patient IDs for EMR/EHR.

- Ultra-thin 1.5-in-high (3.8-cm) platform — slimmest in the industry
- 0.9-in-high (2.3-cm) LCD weight digits for clear readouts
- USB-B port and HL7 IEEE 11073 compliance
- Built-in Wi-Fi and Bluetooth for wireless EMR/EHR connectivity
- Automatically displays full measurements within 1 second of stepping on the platform



Wi-Fi Bluetooth



World-Renowned Weigh Beam Eye-Level Mechanical Physician Scales



The dual-reading height rod measures in inches or centimeters at a glance.

- Available in 13 models with lb, lb/kg, or kg configurations
- Optional accessories: handposts, height rods, and wheels
- Easy-to-clean removable platform cover
- 59-in (150-cm) overall height
- No power required — durable, long-lasting mechanical design



450 lb x 4 oz, 200 kg x 100 g,
or 440 lb x 4 oz & 200 kg x 100 g

Capacity

10.75 in x 14.75 in /
27.3 cm x 37.5 cm

Platform Size

(on select models)
24–84 in / 60–212 cm

Height Rod Range

Dual-reading, die-cast
weighbeam is easy to read
from either side, front or back

Other Features

Mild steel construction with
electrostatic powder-paint finish



 **SINCE
1900**

DETECTO

360-Degree-Rotating Rescue Carts for Mobile Clinical Storage

RESCUE
medical carts



DETECTO

Built with the busy clinical user in mind, Rescue Series carts deliver dependable performance in fast-paced clinical environments, pairing hygienic surfaces with secure, centralized drawer locking.

- New ergonomic handles and push-to-release countertop for comfort and added workspace
- Durable powder-coated steel body with easy-to-clean ABS countertop for hygienic wipe-down disinfection
- 12 standard models available in red, blue, yellow, or white with 5 or 6 drawers; pediatric model includes 9 multi-colored drawers using the Broselow color-coding method for rapid supply identification
- **Cart Dimensions:** 30.5 in x 21.3 in x 42 in / 77.5 cm x 54.1 cm x 106.7 cm
- Central keyed or EMG breakaway lock secures all drawers simultaneously
- 5-in locking casters (2 parking, 1 directional steering), protective ABS corner bumpers, and 360° swivel ensure safe, stable movement through tight spaces
- Auto-return drawer rails and optional accessories for customization or fully-loaded configurations



Technologically Advanced MobileCare Carts with Proprietary Onboard & PC Software

MOBILECARE 
medical carts

MobileCare electronic carts come with a touchscreen LCD for entering pincodes and some models include touchless RFID card badge scanning.

- The most advanced security level features available in the medical cart market
- Advanced customization with nearly unlimited configuration options: colors, handrails, drawer widths (16.5 in/42 cm, 23 in/58 cm, and 29 in/73.5 cm), and a wide range of accessories
- Complementary centralized software manages all employee access and security settings via secure Wi-Fi/Bluetooth for instant updates
- Electronic models include detailed employee audit tracking to help mitigate narcotics theft
- Standard 5-in (12.7-cm) anti-static locking casters and full wrap-around bumper trim for safe, nimble mobility in busy care settings
- 3 secure lock options: electronic (RFID with 3 frequencies + full-color touchscreen pincode), quick-release, or keyed
- Crash cart automatically alerts biomed when used, displays ready/used status via low-power e-paper screen, manages power intelligently, and tracks accessed drawers for targeted restocking



DETECTO



Purpose-Built Pediatric Scale Cart for Streamlined Infant Weigh-Ins

PEDIATRIC
scale cart



DETECTO

Custom-fit design ensures stable, reliable placement of the DETECTO 6745 pediatric scale.

- Flat surface designed to fit most baby scales, including the DETECTO 6745 scale
- Optimized, ergonomic working height reduces caregiver strain during infant weigh-ins and length checks
- Integrated enclosed storage with fixed shelf keeps diapers, wipes, bandages, and other essentials organized and within reach
- Protective edge lip enhances safety and stability by keeping the scale and supplies securely in place during use
- Quiet, hospital-grade locking casters provide smooth, secure mobility throughout patient areas
- Compact, space-saving design combines scale positioning and supply storage into one efficient mobile cart
- Quality aluminum cart body with powder-coat paint and ABS countertop and bumpers
- **Cart Dimensions:** 33.07 in x 21.63 in x 30.94 / 84.01 cm x 54.94 cm x 78.58 cm
- 4-in-diameter (10.16-cm) wheels with 2 parking locks and 1 directional steering lock



8500 Portable Stretcher Scales for Stretcher, Gurney, & Wheelchair Patients




*Model 8550 includes
easily-collapsible column
for weight indicator.*

*Lightweight for easy lifting and
rolling by nursing staff; slim profile
fits through nearly any doorway.*

- Ideal for emergency admissions when rapid, accurate patient weights are critical
- Suited for ER, dialysis, long-term care, and nursing facilities
- Spacious platform accommodates all common stretcher types
- Single-pass weighing allows entry of stretcher or wheelchair tare weight via the indicator's numeric keypad

Wi-Fi Bluetooth™

Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



1,000 lb x 0.2 lb /
450 kg x 0.1 kg

Capacity

60 in x 32.5 in x 2.2 in /
152.4 cm x 82.6 cm x 5.6 cm

Platform Size

6 C batteries (not included)
or AC adapter (included)

Power

Other Features

Two integrated ramps
included for easy access

Column-mounted indicator
or remote display available

Low-profile platform
with optional guide rails

OneWeigh units locking
into pounds or kilograms



DETECTO®

7150-AC ADA-Compliant Wheelchair Scale with Extended Platform & Side Guards



Bluetooth™

Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems

Capacity

1,000 lb x 0.2 lb /
450 kg x 0.1 kg

Platform Size

40 in x 36 in x 2.2 in /
101.6 cm x 91 cm x 5.6 cm

Ramp Size

36 in x 17.4 in (each) /
91.4 cm x 44.15 cm

Power

6 C batteries (not included)
or AC adapter (included)

Other Features



BMI calculation

OneWeigh units locking into
pounds or kilograms

Easily-collapsible column
with quick release lever
and safety pin

Ships 90% assembled

- Meets ADA accessibility compliance standards with two-way, 1:8 shallow-angle ramps, wider platform, and side guards that help prevent wheelchairs from rolling off the scale
- Vertical fold-up design and integral wheels enable easy transport and space-saving storage
- Up to 99 optional tares allow storage of wheelchair weights for faster, more efficient workflows



DETECTO

By August 9, 2026, under the ADA's Medical Diagnostic Equipment (MDE) accessibility requirements, certain healthcare facilities must have at least 10% accessible medical equipment, but no fewer than one accessible weight scale (20% for mobility rehabilitation facilities). Facilities should consult applicable regulations to determine compliance requirements.



6550 Portable Wheelchair Scale with Two-Way Built-In Ramps



With only one hand you can quickly slide the release lever and fold down the 6550's column down for easy transport or vertical storage.

- Vertical fold-up design enables space-saving storage
- Up to 99 optional tares allow storage of wheelchair weights for faster, more efficient workflows
- Convenient pushbutton or keypad tare quickly subtracts wheelchair weight for accurate patient readings



1,000 lb x 0.2 lb /
450 kg x 0.1 kg

Capacity

32 in x 36 in x 2.2 in /
81 cm x 91 cm x 5.6 cm

Platform Size

6 C batteries (not included)
or AC adapter (included)

Power

BMI calculation

Other Features

Two-way, low-profile platform for
easy patient access from either side

Easily-collapsible column with quick
release lever and safety pin

OneWeigh units locking into
pounds or kilograms



Bluetooth™ Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



6560 Multi-Use Portable Wheelchair Scale with Curved Handrail



Bluetooth™

Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



Fold-up seat included with 6570 model for seated, standing, and wheelchair weighing.

Capacity

1,000 lb x 0.2 lb /
450 kg x 0.1 kg

Platform Size

32 in x 36 in x 2.2 in /
81 cm x 91 cm x 5.6 cm

Power

6 C batteries (not included)
or AC adapter (included)

Other Features

Elegant, curved handrail
for patient stability while
weighing standing

May be used for standing
or wheelchair weighing
(or seated with model 6570)

Column folds down for easy
transport or vertical storage

Integral wheels ease transport

OneWeigh units locking
into pounds or kilograms

BMI calculation



- Dual-access ramps extend fully to the floor, allowing maximum convenience for loading and unloading patients from either side
- Unique ramp design requires half the effort to roll a wheelchair-bound patient onto the scale compared to competitor models
- Low-profile platform enhances patient accessibility and ease of use
- Convenient pushbutton or keypad tare quickly subtracts wheelchair weight for efficient workflows



DETECTO®



BRW1000 Portable Bariatric Wheelchair Scale with Spacious Low-Profile Platform

The MedVue® multifunctional indicator with alphanumeric keypad, RS232 serial, and USB ports displays height, weight, and BMI onscreen simultaneously. It's capable of 14-digit patient ID entries and 99 stored wheelchair tare weights.



Bluetooth™



Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



- Accommodates bariatric wheelchairs up to 40 in (101.6 cm) wide for maximum accessibility
- Chair tare function allows quick entry of wheelchair weight for efficient, comfortable patient weighing
- Heavy-duty platform with non-skid rubber mat supports manual chairs, power chairs, dialysis patients, and walk-ons
- 4 stainless-steel load cells with self-aligning feet automatically level the scale after transport for accurate measurements



1,000 lb x 0.2 lb /
450 kg x 0.1 kg

Capacity

32 in x 40 in x 2.2 in /
81 cm x 101 cm x 5.6 cm

Platform Size

6 C batteries (not included)
or AC adapter (included)

Power

BMI calculation

Other Features

Two-way, easy access,
integral ramps

Integral wheels and
guide handles

OneWeigh units locking
into pounds or kilograms

Flexible coiled cable
remote display extends up
to 14 ft (427 cm) from the scale



DETECTO



6400 Portable Digital Wheelchair Scale with Two-Way Easy-Access Ramps



Bluetooth™

Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



Capacity

1,000 lb x 0.2 lb /
450 kg x 0.1 kg

Platform Size

32 in x 36 in x 2.2 in /
81 cm x 91 cm x 5.6 cm

Power

6 AA batteries (not included)
or AC adapter (included)

Other Features

Flexible coiled cable remote display extends up to 14 ft (4.2 m)

Integral dual wheel and guide handle for transport

Low-profile platform for easy patient access

Vertical mode for transport and storage

OneWeigh units locking into pounds or kilograms

BMI calculation

Model 750 indicator with easy-to-use, basic weighing functions, RS232 serial port, and big bold LCD digits.



- Built-in entry and exit ramps support all wheelchair weighing applications
- Low platform height enhances patient accessibility and comfort
- Platform and ramps feature black rubber matting for secure, non-slip traction
- Versatile indicator bracket allows wall mounting for stationary use or tabletop placement for mobile operation


DETECTO



7550 Connected Wall-Mount, Fold-Up Stationary Wheelchair Scale



Integrated full-color, touchscreen LCD with wheelchair tare weight mode to remove wheelchair weights and factor patient weight.



1,000 lb x 0.2 lb /
450 kg x 0.1 kg

Capacity

36 in x 32 in x 1.5 in /
91.4 cm x 81.3 cm x 3.8 cm

**Platform Size
(Flat Portion)**

6 C batteries (not included)
or AC adapter (included)

Power

BMI calculation

Other Features

- Wall-mount design installs easily in high-traffic areas and stores neatly against the wall
- Fold-up platform saves valuable floor space and keeps pathways clear when not in use
- Gas-spring mechanism gently lowers the platform into the weighing position
- Slim, non-weighing mode blends in for an efficient workflow

Slim 7-in (17.8-cm) profile minimizes space usage in hallways, clinical offices, and patient rooms

Preserves patient dignity by allowing weighing in open common areas

Integrates into any healthcare environment with unique fold-up storage

OneWeigh units locking into pounds or kilograms



DETECTO



6855 Waist-High Digital Bariatric Scales with Wraparound Handrails



Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



Capacity

600 lb x 0.2 lb /
270 kg x 0.1 kg

Platform Size

18 in x 14 in x 2.4 in /
46 cm x 36 cm x 6.1 cm

Power

6 AA batteries (not included)
or AC adapter (included)

Other Features

Stainless steel handrails
with padded grips for
patient assurance

Integral wheels for quick
and easy transport

Digital and mechanical
height rods available

OneWeigh units locking
into pounds or kilograms

BMI calculation

Model 6854DHR
with digital height
rod that measures
up to 79 in (200 cm).

Model 6855MHR with
mechanical height rod
that measures up to
84 in (212 cm).

- Designed to safely and easily weigh obese or unsteady patients
- Low-profile, spacious platform with non-skid surface enhances patient comfort and stability
- Weight-lock function allows clinicians to assist patients off the scale before recording results

DETECTO



6856 Waist-High Digital Bariatric Scales with Wraparound Handrails

Model 6857DHR with digital height rod that measures up to 79 in (200 cm).

Model 6856MHR with mechanical height rod that measures up to 84 in (212 cm).

- Balanced design with integral wheels for easy movement between rooms
- Low-profile, spacious square platform, perfect for measuring obese or unsteady patients
- Weight-lock function allows clinicians to assist patients off the scale before recording results



Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



1,000 lb x 0.2 lb /
450 kg x 0.1 kg

Capacity

24 in x 24 in x 2.4 in /
61 cm x 61 cm x 6.1 cm

Platform Size

6 AA batteries (not included)
or AC adapter (included)

Power

Stainless steel handrails
with padded grips for
patient assurance

Other Features

Integral wheels for quick
and easy transport

Digital and mechanical
height rods available

OneWeigh units locking
into pounds or kilograms

BMI calculation



DETECTO

6868 Bariatric Flip Seat Digital Scale with Wraparound Handrails



Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



Model 750 indicator with easy-to-use, basic weighing functions, RS232 serial port, and big bold LCD digits.

Capacity	1,000 lb x 0.2 lb / 450 kg x 0.1 kg
Platform Size	32 in x 33.7 in x 2.4 in / 81.3 cm x 85.7 cm x 6 cm
Power	6 AA batteries (not included) or AC adapter (included)
Other Features	Unique flip-up seat for standing or seated weighing Integral wheels and balanced design for mobility Padded stainless steel wrap-around handrails for secure, comfortable patient support Automatic shutoff/sleep modes OneWeigh units locking into pounds or kilograms BMI calculation



Rubber mat provides traction and stability.



- Versatile 3-in-1 design supports standing, seated, and bariatric weighing
- Portable design with integral wheels and battery power for easy transport
- Lock/Release key holds weight values so clinicians can assist before recording
- Slim profile fits through doorways up to 32 in (81 cm) wide

DETECTO



6800 Portable Bariatric Floor Scales with Wall/Desk-Mounted Indicator



Bluetooth™ Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



1,000 lb x 0.2 lb /
450 kg x 0.1 kg

**Capacity
(6800)**

300 kg x 0.1 kg
or 350 kg x 0.2 kg

**Capacity
(6800KGEU)**

24 in x 24 in / 61 cm x 61 cm

Platform Size

6 AA batteries (not included)
or AC adapter (included)

Power

BMI calculation

Other Features

Integral carrying handle

Low-profile design

OneWeigh units locking
into pounds or kilograms



- Ideal for clinical weighing applications, with a low-profile platform is designed to enhance patient comfort
- Includes universal bracket and 10-ft (3-m) indicator cable for flexible indicator mounting
- Skid-proof, textured rubber surface for improved patient stability
- Lightweight design for easy portability



DETECTO

6475 Portable Digital Chair Scales with Lift-Away Armrests



Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



Both armrests lift for easy lateral access from a wheelchair or bed, and the fold-up footrest provides added support during weighing while retracting neatly when not needed.

Capacity

400 lb x 0.2 lb /
180 kg x 0.1 kg

Seat Dimensions

16.6 in x 16.7 in x 14.3 in /
42.2 cm x 42.4 cm x 36.3 cm

Power

6 AA batteries (not included) or AC adapter (included)

Other Features

Lift-away armrests and fold-up footrest

Plastic-molded seat for patient comfort

Integral wheels and guide handle ease transport

Sturdy steel construction

OneWeigh units locking into pounds or kilograms



Analog model available (model 475).



- Arrive fully assembled for immediate use
- Ideal for disabled, dialysis, and elderly patients who need seated weighing
- Comfortable, stress-free measurement of weight and BMI

DETECTO®



6880 Rolling Digital Chair Scale with Smart-Phone-Style Indicator



Smartphone-style indicator with bright 0.75-in-high (1.9-cm) LCD, simple 6-button keypad, battery level indicator, and temporary weight-lock function for clear viewing and recording.



Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



500 lb x 0.2 lb /
250 kg x 0.1 kg

Capacity

16.6 in x 16.7 in x 14.3 in /
42.2 cm x 42.4 cm x 36.3 cm

Seat Dimensions

12 AA batteries (not included)
or AC adapter (included)

Power

Other Features

Contoured, sloped backrest
promotes a relaxed
patient position

Padded armrests enhance
comfort and fold away
for easy access

Dual integral wheels
and guide handles for
easy maneuverability

Stylish black paint color
that blends into any
healthcare environment

OneWeigh units locking
into pounds or kilograms

- Designed for patients with limited mobility
- Ideal for nursing homes, hospitals, geriatrics, and orthopedic departments
- Ergonomic seat provides comfortable, stable weighing
- Intuitive interface streamlines operation for healthcare staff
- Retractable armrests and footrest for simplified patient transfers



DETECTO

6745 Digital Baby Scales with 1-Gram Divisions for Milk Intake & Wet Diaper Weighing



Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems

Features 5 stored memory weights plus Milk and Diaper weighing functions for tracking infant fluid output.

Capacity

0–10 lb x 0.002 lb (0.02 oz) /
10–30 lb x 0.005 lb (0.05 oz);
0–5 kg x 0.001 kg
(0–5,000 g x 1 g) /
5–15 kg x 0.002 kg
(5,000–15,000 g x 2 g)

Platform Size (Flat Portion)

22.5 in x 9.1 in /
57.1 cm x 23.1 cm

Power

6 C batteries (not included)
or AC adapter (included)



Other Features

OneWeigh units locking
into pounds or kilograms

Full-color touchscreen LCD

Integral analog measuring
tape up to 22 in (56 cm)

RS232 for printer output
and USB-B for PC connectivity



- High-precision design for NICU, nursery, delivery room, and pediatric use
- NTEP certification ensures verified accuracy and helps prevent dosing errors
- Lift-off clear polystyrene tray allows fast, hygienic cleaning
- Built-in carrying handles for easy transport
- Wide, stable base enhances infant safety during weighing
- Colorful, child-friendly graphics make weighing fun and attractive for young patients

DETECTO



8432-CH Digital Baby Scales with Inclined Plastic Chair-Like Seat

The weight indicator is conveniently positioned above the seat for clear, easy-to-read measurements while the baby is seated.



- Easy-to-clean inclined plastic seat with a durable enameled steel base keeps infants and toddlers comfortable and secure
- Seashell-shaped seat design makes weighing fun and relaxing for pediatric patients
- Dual pound/kilogram readouts for versatile clinical use
- Heavy-duty steel understructure provides stability and long-lasting performance



44 lb x 0.05 lb /
20 kg x 0.02 kg

Capacity

26 in / 66 cm

Scale Height

6 AA batteries (not included)
or AC adapter (included)

Power

OneWeigh units locking
into pounds or kilograms
for safe dosages

Other Features

Solid steel base for steady
assurance while weighing

Pushbutton tare removes
blanket weights

DETECTO



8440 Digital Baby Scale with Collapsible Tray for Standing Toddler Use & Easy Transport



Collapsible tray allows compact storage, easy transport, and quick conversion for standing toddler weighing.



Capacity	44 lb x 0.5 oz / 20 kg x 10 g
Tray Dimensions	20.6 in x 10.5 in x 4 in / 52.4 cm x 26.7 cm x 10.2 cm
Power	1 9V battery (not included) and 1 3V CR2032 lithium back-up battery (included) or AC adapter (included)
Other Features	Storage of previous weight memory recall Calculation of change in weight Large, easy-to-read display Removable weighing tray



Optional carrying case available (model 8440-CASE).



- Accurate, safe weighing for infants and toddlers—ideal for home health and clinics
- 1-in-high (2.5-cm) LCD displays pound/kilogram readings and includes a built-in clock
- Weight Difference key compares current weight to a stored value with time/date stamp

DETECTO



MB130 Digital Baby Scale with Extra-Large Removable Tray for Standing Toddler Use



Optional digital length measuring rod (model DLM) mounts easily to the scale for a fully integrated, all-digital measuring station.



Optional carrying case with adjustable strap available (model MB-CASE).

- Dual graduations for neonatal-to-pediatric use
- Removable tray converts the unit into a low-profile step-on toddler scale
- 1-in (2.5-cm), 5-digit LCD ensures clear, highly visible readouts
- Comfortable, contoured weighing tray cradles babies securely
- Portable, lightweight, and battery-powered for mobile medical use



0–22 lb x 0.01 lb /
22–44 lb x 0.02 lb;
0–10 kg x 5 g /
10–20 kg x 10 g

Capacity

22 in x 11.3 in x 3.7 in /
55.9 cm x 28.7 cm x 9.4 cm

Tray Dimensions

1 9V battery (not included)
or AC adapter (included)

Power

Integral analog measuring tape

Other Features

Easy to use with large readouts

Auto shut-off after 2 minutes
to conserve battery life

Breast Milk Intake function
compares baby weight
before and after feeding



DETECTO



DLM Digital Length Measuring Rod for Infants



Convenient unit switching: choose from ft/in, in, or cm.




Length Range	13.75–31.5 in / 35–80 cm
Graduations	0.0625 in / 1 mm
Dimensions	11.4 in x 2.7 in x 24.4 in / 29 cm x 7 cm x 62 cm
Power	2 1.5V AAA batteries (included)
Other Features	Easy-to-read LCD display

Hold key retains measurements for active babies

Switchable units

Smooth, contoured paddles are gentle to infant touch

Easily mounts to model MB130.



- May be used stand-alone on an exam table or attached to the MB130 baby scale
- Lightweight, battery-powered design for easy room-to-room transport
- Head and foot positioners ensure accurate infant and toddler length measurements
- Bold LCD with convenient unit switching for clear, versatile readouts
- Rubber feet provide stable, slip-resistant positioning during use

DETECTO



PHR Free-Standing Mechanical Portable Height Rod




PHR works in conjunction with DR400C for simultaneous height/weight measurement.



Optional carrying case with shoulder strap available (model PHR-CASE); fits the PHR and DR400C.

- May be used stand-alone or in conjunction with the DR400C scale
- Simple, quick setup with no wall fastening required
- Can be used free-standing or placed against a wall; large floor plate for stability
- Includes two wall spacers for added stability when positioned against a wall
- Reversible head piece allows pairing with DR400C scale for simultaneous height/weight measurement; head piece flips upside down when used without scale



4.5–81 in / 11.5–205 cm

0.125 in / 0.1 cm

13 in x 19.8 in x 83.5 in /
33 cm x 50.3 cm x 212.1 cm

From pediatric to
adult measurement

Height rod disassembles
easily into 4 sections and
stores conveniently with
head piece and base


Durable and lightweight
all-plastic construction

Height Range

Graduations

Dimensions

Other Features



DETECTO



SONARIS Touchless Sonar Height Rod with State-of-the-Art Triangulation Technology



Height Range	0–9 ft / 30–275 cm
Graduations	0.2 in / 0.5 cm
Dimensions	5.9 in x 3.5 in x 47.4 in / 14.9 cm x 8.9 cm x 120.3 cm
Power	6 AA batteries (not included) or AC adapter (not included)
Other Features	4 units of measurement: ft/in, in, cm, and mm

Instant digital height readouts
without pushing a single button

Auto or manual calibration
to current position for
full portability



DETECTO



*Easily toggle between 4
units of measurement:
ft/in, in, cm, or mm.*



- Ideal for school health and medical offices
- Compact readout with 0.75-in (1.9-cm) backlit LCD for clear visibility; desk or wall mountable
- Can be used stand-alone on a tabletop (or any flat surface) for true mobility



Solace Series Flush-Mounted Platform Scales for Wheelchair & Stretcher Weighing



Multi-color touchscreen indicator with built-in printer displays weight, height, and BMI simultaneously with a clear 0.9-in-high (2.3-cm) LCD; includes auto weight-lock for easy recording.



1,000 lb x 0.2 lb /
450 kg x 0.1 kg

Capacity

7 standard sizes:

36 in x 32 in / 91 cm x 81 cm,
36 in x 36 in / 91 cm x 91 cm,
48 in x 36 in / 122 cm x 91 cm,
48 in x 48 in / 122 cm x 122 cm,
72 in x 36 in / 183 cm x 91 cm,
72 in x 48 in / 183 cm x 122 cm,
90 in x 48 in / 229 cm x 122 cm

Platform Size

Corrosion-resistant
stainless steel scale deck

Other Features

Pit frame, trim ring, and
wall-mount indicator
combo included

Printer and handrail
are optional

OneWeigh units locking
into pounds or kilograms

- 28 standard models to choose from (custom sizes/configurations available)
- Ideal for dialysis monitoring and bariatric weighing
- In-floor, stainless-steel platform installs flush with the floor to remove tripping hazards and improve accessibility
- Accommodates walk-on or wheel-on use; 6-ft (183-cm) platform fits stretchers and beds — no dual-scale setup needed
- 110 VAC digital indicator provides weight readouts in pounds or kilograms



APEX-RI Portable Digital Scales with Remote Indicator



High-tech, smart-phone-style remote indicator with 6 easy-use buttons, bold 0.75-in-high (1.9-cm) clinical-blue LCD readouts, and a 6-ft (183-cm) cable for flexible viewing.



Capacity 600 lb x 0.1 lb / 270 kg x 0.05 kg

Platform Size 17 in x 17 in x 2.75 in /
43 cm x 43 cm x 7 cm

Power 12 AA batteries (not included)
or AC adapter (included)

Other Features Mother/Baby function subtracts
adult weight to capture infant/
toddler weight while held

Wall-mount or desktop
mounting bracket included

Integral carrying handle

OneWeigh units locking
into pounds or kilograms

BMI calculation



- Extra-wide, low-profile platform with textured surface for weighing traction and user safety
- 1 RS232 serial port and 1 micro USB-B port; HL7 IEEE 11073 compliant
- Ready to use out of the box with no assembly needed



DETECTO



Wi-Fi/Bluetooth included on -C
models for wireless connectivity
with EMR/EHR systems



DR660 Portable Digital Bariatric Scale with Rechargeable Built-In Battery Pack

Integral carrying handle and lightweight design for easy transport; coiled scale-to-indicator cable extends 34–126 in (86–320 cm).



Adjustable tilt display, bold 1-in (2.5-cm) LCD readouts, HOLD/TARE keys, and lb/kg units locking for easy, accurate weighing.



660 lb x 0.5 lb /
300 kg x 0.2 kg

Capacity

21.5 in x 17.2 in x 1.97 in /
56 cm x 44 cm x 5 cm

Platform Size

RS232 serial output

Other Features

Integral carrying handle

Rubber non-skid mat

Remote display may be mounted on a wall or placed on a tabletop

- Ideal for home healthcare and bariatric weighing
- Built-in rechargeable battery provides up to 8 hours of continuous use; AC adapter included

- Adjustable feet for stability on various weighing surfaces
- Broad platform with low-profile step for easy access
- Optional automatic shutoff helps conserve battery life



DETECTO

GP-600-MV1 Portable Digital Low-Profile Scales for Athletic Wrestling or Dietary Clinic Weigh-Ins



Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



Vivid 3.5 in x 3 in (8.9 cm x 7.6 cm) backlit LCD with nearly 1-in-high (2.2-cm) readouts for easy visibility.



Capacity	600 lb x 0.2 lb / 270 kg x 0.1 kg
Platform Size	18 in x 14 in x 2 in / 46 cm x 36 cm x 6 cm
Power	6 C batteries (not included) or AC adapter (included)
Other Features	Bubble level on platform Integral carrying handle OneWeigh units locking into pounds or kilograms BMI calculation



- Spacious, easy-to-clean platform with black non-skid mat and bubble level for safety assurance
- Lock/Release key holds weight after stepping off the scale
- USB and RS232 connectivity standard

DETECTO



DR400C & DR550C Portable Digital Scales, Compact for Space-Saving Efficiency



Raised, tactile keys with audible feedback and a bright backlit LCD with bold 1-in-high (2.5-cm) digits for clear visibility in any lighting.



- Arrive fully assembled and ready to use out of the box
- Lightweight, low-profile platform for easy portability
- Ideal for mobile clinics, home healthcare, physician offices, and general-purpose/home weighing
- Indicator includes a versatile drop-in mounting bracket with coiled cable for floor, tabletop, or wall-mounted viewing
- Slip-resistant textured rubber mat and non-skid feet for secure traction and stable weighing
- Optional padded carrying case available



400 lb x 0.5 lb /
180 kg x 0.2 kg

**Capacity
(DR400C)**

550 lb x 0.2 lb /
250 kg x 0.1 kg

**Capacity
(DR550C)**

12.5 in x 11.9 in x 2.2 in /
31.8 cm x 30.3 cm x 5.8 cm

Platform Size

4 AA batteries (not included)
or AC adapter (included)

Power

Backlit blue LCD with
bold 1-in-high readouts

Other Features

Coiled cable from base for
versatile indicator positioning



DETECTO

SlimTALK Digital Talking Bathroom Scales with Bilingual Spoken Weight Announcements



The SlimTALK XL is also available with a higher 550 lb / 250 kg capacity and a ribbed black mat.

Capacity	400 lb x 0.1 lb / 180 kg x 0.1 kg
Dimensions	14.3 in x 12.3 in x 1.4 in / 36.3 cm x 31.2 cm x 3.6 cm
Power	3 AAA batteries (included)
Other Features	Speaks 2 languages: English and Spanish Auto-on/off function Durable plastic platform with textured surface for weighing traction and user safety



The step-on/off function enables weighing without stooping to turn on the scale and helps conserve battery life.

- Ideal for elderly, obese, blind, and telemedicine use
- Easy to use for the visually impaired with clear, pleasant spoken announcements of weight (English/Spanish or Off)
- Extra-wide, low-profile platform provides easy step-on access for users of all sizes
- Large 1.5-in-high (3.8-cm) LCD digits for easy viewing

DETECTO



SlimPRO Digital Talking Bathroom Scales with Bilingual Spoken Weight Announcements



Mother/Baby mode subtracts adult weight to capture infant or toddler weight while held.



600 lb x 0.1 lb /
272 kg x 0.1 kg

Capacity

17 in x 14.5 in x 1.4 in /
43 cm x 37 cm x 3 cm

Dimensions

4 AAA batteries (included)

Power

Backlit blue LCD

Other Features

Integral carrying handle

Mother/Baby weighing mode

Non-skid ribbed rubber mat
for weighing traction
and user safety

Speaks 2 languages:
English and Spanish

Auto-on/off function

BMI calculation



- Extra-wide, low-profile platform provides easy step-on access for users of all sizes
- Easy to use for the visually impaired with clear, pleasant spoken announcements of weight (English/Spanish or Off)
- Large 1.5-in-high (3.8-cm) digits on a backlit blue LCD for easy reading
- Stylish gray color blends into any environment

DETECTO



WPS12UT Submersible Washdown Digital Stainless Steel Scale with Utility Bowl



*Submersible underwater
for hygienic cleaning.*



Capacity

12 lb x 0.1 oz / 12 lb x 0.005 lb /
12 lb x 0.125 oz / 192 oz x 0.05 oz /
192 oz x 0.125 oz / 5,500 g x 1 g

Dimensions (Scale and Bowl)

13.6 in x 13.6 in x 6.9 in /
34.5 cm x 34.5 cm x 17.5 cm


Other Features

Accurately measures fluid
or blood loss during clinical
and procedural use

Excellent for weighing diaper
fluids or blood absorbed by
chux in labor and care settings

Tare function removes dry
diaper or chux weight for
precise fluid measurement

NSF approved



*Stainless steel diaper tray
available with model WPS12DT.*



- IP67 rated, fully submersible for easy washdown cleaning
- Ideal for weighing wet diapers, used chux underpads, lap sponges, and organs
- 80-hour rechargeable battery with automatic shut-off for portable use (AC adapter also included)
- 4.6-in (11.7-cm) deep stainless steel utility bowl for sample retention and easy cleanup
- Backlit LCD with large 1.2-in (3-cm) readouts for clear visibility in low light
- Raised, tactile buttons for easy, confident operation



DETECTO



420 Series Digital Precision Balance Scales for Retail, Laboratory, & Medical Weighing



Nearly 1-in-high (2.3-cm) bold digits with blue backlit LCD.



Excluding models 420-100 and 420-2000

- High-precision weighing with accuracy down to 10 mg for demanding applications
- Removable air shield included on select models for enhanced measurement accuracy
- Built-in rechargeable battery for portable, uninterrupted operation
- Ideal for retail applications including herbs, cannabis, specialty coffee and tea, jewelry, and precious metals
- Hospital-ready performance for weighing blood bags, blood components, patient medical waste samples, and pill counting in hospital pharmacies
- Laboratory versatility for weighing chemicals and powders in educational, industrial, and agricultural settings



7.4 in x 5.8 in /
18.8 cm x 14.8 cm

**Tray Size
(420-3000)**

5.7 in x 4.5 in /
14.5 cm x 11.5 cm

**Tray Size
(all other models)**

Accuracy to 10 mg

Other Features

8 different units of measure

Class II legal for trade

Counting mode

RS232 serial port
for printer output

Polycarbonate body
with stainless steel tray

Air shield included with
all models except 420-3000



DETECTO




IB Series In-Bed Scales with Hydraulic Lift Mechanism for Weighing Bed-Bound Patients



Wi-Fi/Bluetooth included on -C models for wireless connectivity with EMR/EHR systems



Hydraulic lift design simplifies weighing of bedridden, non-ambulatory patients.



Capacity (IBFL500)

500 lb x 0.1 lb / 225 kg x 0.1 kg

Capacity (IB800)

800 lb x 0.2 lb / 350 kg x 0.1 kg

Power


6 AA batteries (not included) or AC adapter (included)

Other Features

Stretcher carrying strap and "easy stow" hoops

Full-swivel ball-bearing casters with rear locking pair for secure weighing

- Flame-retardant, antibacterial 6-ft-long (1.8-m) stretchers placed beneath the patient and hydraulically lifted from the bed for weighing
- Sturdy steel construction with locking casters and an adjustable weighing arm
- ZERO key subtracts bedding weight for accurate measurements
- **IB800**: adjustable base; fits through 32-in (81-cm) doorways
- **IBFL500**: fixed base; fits through 36-in (91-cm) doorways
- Not intended for patient transport



DETECTO®



PL Series Digital Patient Lift Scales with Portable Lift Mechanisms for Suspended Weighing



400 lb x 0.2 lb / 180 kg x 0.1 kg

**Capacity
(PL400)**

600 lb x 0.2 lb / 270 kg x 0.1 kg

**Capacity
(PL600)**

4.175 in x 2.175 in x 5 in /
10.6 cm x 5.52 cm x 12.7 cm

Dimensions

1 9V battery (not included)

Power

Suspended weighing

Other Features

Lock-in weight option

0.75-in-high (1.9-cm) LCD readouts

- Suspended weighing for patients lifted by portable patient lift mechanisms
- Simple 3-key membrane keypad (ON/ZERO, OFF, lb/kg) for easy operation
- Includes 2 spring O-rings (top and bottom)

- Lock-in weight feature stores the last reading when powered off
- Battery powered with automatic shut-off to conserve battery life



DETECTO

Step-On Cans: Ruggedly-Built, Flame-Retardant Metal Waste Receptacles Available in 3 Finishes



Model PR Series

Model P Series

Model C Series

Model P-32R.

- Rugged, flame-retardant construction ideal for hospitals, doctors' offices, patient rooms, and laboratory environments
- Heavy-duty foot pedal and linkage with self-closing lid for hands-free, sanitary operation
- Unitized metal construction with integrated bag-securing mechanism
- Available in 5 sizes (16 qt/15 L, 24 qt/23 L, 32 qt/30 L, 48 qt/45 L, or 100 qt/95 L) and 3 finishes (stainless steel; baked epoxy painted red or white)
- Biohazard and infection-waste adhesive labels included



DETECTO



Excluding models C-100,
P-100, and P-100R



Waste Mate Commercial-Grade Plastic Waste Receptacles in 4 Sizes & 2 Standard Colors



Foot pedal flips open lid for quick trash disposal.

- Designed for heavy-duty use in any fast-paced environment
- Hands-free foot pedal opens the lid for no-touch disposal; soft-close lid lowers quietly into place
- Integrated internal clamps and quick-replacement hinges secure trash bags and allow fast bag changes
- High-grade plastic construction resists rust, scratches, and denting and is easy to clean
- Reinforced rim and steel wire outriggers (front and back) enhance stability and reduce tipping during pedal operation
- Slim profile (12.8 in/32.4 cm or 15.5 in/39.3 cm) ideal for tight spaces
- Ergonomic upper lip for comfortable handling and easy transport
- Available in 4 sizes (8–24 gallons) and 2 colors (red and beige)
- Limited-angle lid opening prevents a full 90° swing, protecting nearby walls and furniture from damage




DETECTO®



T3530 Series Hanging Dial Scales with Removable Stainless Steel Weighing Bowl



Capacity	32 lb x 1 oz / 15 kg x 25 g
Dial Size	14.5-in / 37-cm diameter
Other Features	Readable from both sides, front and back Stainless steel scale and bowl construction



Model T3530 with hanging scale stand (model HS-STAND), sold separately

Use Model T3530 for avoirdupois weighing.



Use Model T3530KG for metric weighing.



- Deep removable stainless steel bowl (13.5-in / 34-cm diameter) with 3 drainage holes, ideal for foodservice and medical weighing, including autopsy use
- Extra-large, easy-to-read dial with protective glass plates for durability
- Dual-sided dial visibility for convenient reading from either side
- Available in avoirdupois or metric models
- Thermolastic springs ensure long-term accuracy despite extreme temperature variations
- Air dash pot shock absorption protects internal mechanisms and extends scale life

DETECTO



GH Series Universal Glove Box Holders in Mild Steel & Stainless Steel



- Space-saving universal design conserves counter and shelf space
- Holds 1, 2, or 3 glove boxes depending on model
- Open-face construction allows easy identification of glove size and type
- Available in 3 sizes and 2 finishes (white baked epoxy or stainless steel)
- Backed by DETECTO lifetime guarantee



6.5 in x 4.25 in x 10.5 in /
17 cm x 11 cm x 27 cm
(models GH1 and GH1SS)

12 in x 4 in x 10 in /
30 cm x 10 cm x 25 cm
(models GH2 and GH2SS)

18 in x 4 in x 10 in /
46 cm x 10 cm x 25 cm
(models GH3 and GH3SS)

Dimensions

Compatible with
most glove box types

Other Features

Open face design

Installs in seconds with
mounting hardware included

Private labeling available

Sturdy steel construction







Weigh. Measure. Store. Dispose.

NOTES



102 East Daugherty Street, Webb City, MO 64870 USA
Toll-Free: (800) 641-2008 | Email: detecto@cardet.com
www.detecto.com

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.